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

TIMBRO E FIRMA DEL CLIENTE

Rozzi

Chapter 1 Installation and calibration

This chapter contains all information necessary for installation of a new XW8600 **VolumeShare 4 / Seno Advantage 2.2** product, as well as upgrading an existing XW8400 workstation of the Installed Base:

1 Seno advantage workstation configuration form

Complete this form at the end of installation.
Leave it attached to the site's Service Manual. Upgrade it for each new option or software release.

HARDWARE TYPE:	XW8600	XW8400
SERIAL NUMBER:	C7C03409YL	
HOSTNAME:	RSW	
IP ADDRESS:	
NETMASK:	
DEFAULT ROUTER IP ADDRESS:	
DIRECT-CONNECT IP ADDRESS:	
DIRECT-CONNECT NETMASK	
ROOT PASSWORD:	
SDC PASSWORD:	ADW4.5	
SDC_ADMIN PASSWORD:	
LICENSE ID (licenseld command):	e951e12d 3a6dc4ce	
HOSPITAL NAME & ADDRESS:	OSPEDALE S. MARTINO VIA ROCKEFELLER ORISTANO	
CUSTOMER NAME & EMAIL:	
SYSTEM ORDER (GON ; SO nr):	
SYSTEM ID (SysID):	A5302319	
APPLICATIONS SUPPORT PHONE Nr:	
SALES SUPPORT PHONE Nr:	
USER INTERFACE LANGUAGE:	ITALIAN	
KEYBOARD LANGUAGE:	US	
AW SOFTWARE INITIAL RELEASE:	AW4.5 e951e12d _02.113.ctt-5.X	
AW SOFTWARE UPGRADE RELEASE:	
AW SOFTWARE LICENSE KEY:	e951e12d	
SA SOFTWARE INITIAL RELEASE:	SA22 3.1.1	
SA SOFTWARE UPGRADE RELEASE:	
SA SOFTWARE LICENSE KEY:	232593453A	
PREMIUM VIEW LICENSE KEY:	3077541695	
FLOATING LICENSE KEY:	
VOLUME VIEWER LICENSE KEY:	8380DJSH44ME9PR9 1	
VOLUME VIEWER PET LICENSE KEY:	

Installation and calibration

- AUTOLAUNCH LICENSE KEY:
- AVA XPRESS LICENSE KEY:
- VESSEL IQ XPRESS LICENSE KEY:
- ADVANTAGE CTC PLUS LICENSE KEY:
- ADVANTAGE CTC PRO LICENSE KEY:
- ADVANTAGE CTC PRO 3D LICENSE KEY:
- COLON VCAR LICENSE KEY:
- ADVANTAGE ALA LICENSE KEY:
- LUNG VCAR LICENSE KEY:
- CARDIQ XPRESS 2.0 PLUS LICENSE KEY:
- CARDIQ XPRESS 2.0 PRO LICENSE KEY:
- CARDIQ XPRESS 2.0 ELITE LICENSE KEY:
- CARDIQ FUSION PET LICENSE KEY:
- CARDIQ FUSION SPECT LICENSE KEY:
- CARDEP LICENSE KEY:
- AUTOBONE LICENSE KEY:
- AUTOBONE XPRESS LICENSE KEY:
- In ROOM 3D LICENSE KEY:
- PET VCAR LICENSE KEY:
- CARDIQ FUNCTION XPRESS LICENSE KEY:
- THORACIC VCAR LICENSE KEY:
- VOLUME VIEWER MR LICENSE KEY:
- VOLUME VIEWER BASE LICENSE KEY:
- INTEGRATED REGISTRATION LICENSE KEY:
- DYNAMIC SHUTTLE LICENSE KEY:
- AWCCOW Enabler (enab) 2102600389 (switch) 1754664061
- DENTASCAN LICENSE KEY:
- DENTASCAN PLUS LICENSE KEY:
- CT PERFUSION LICENSE KEY:
- ADVANTAGE 4D LICENSE KEY:
- ADVANTAGE SIM MD LICENSE KEY:
- SMARTSCORE LICENSE KEY:
- CT MR FUSION LICENSE KEY:
- CT PET FUSION LICENSE KEY:
- FUNCTOOL PERFORMANCE LICENSE KEY: 906507389
- FUNCTOOL SPECTROSCOPY LICENSE KEY:
- FUNCTOOL DIFFUSION TENSOR LICENSE KEY:
- FUNCTOOL DIFFUSION TENSOR TRACTOGRAPHY LICENSE KEY:
- FUNCTOOL BRAINSTAT LICENSE KEY:
- BRAINWAVE LICENSE KEY:
- CADSTREAM LICENSE KEY:

MR PASTING LICENSE KEY:

.....

ADVANTAGE PASTE LICENSE KEY:

STENOSIS ANALYSIS LICENSE KEY:

LEFT VENTRUCULE ANALYSIS LICENSE KEY:

.....

REMOTE ACCESS LICENSE KEY:

.....

NETWORK HOST #1: Type:.....

Hostname:.....IP Address:.....Netmask:.....

Protocol:.....Port Number:.....AE Title:.....

Direct Connect IP Address.....Direct Connect Netmask:.....

Query/Retrieve supported:Y.....N

Storage Commitment :

SC IP address:.....SC Port Number:.....SC AE Title:.....

NETWORK HOST #2: Type:.....

Hostname:.....IP Address:.....Netmask:.....

Protocol:.....Port Number:.....AE Title:.....

Direct Connect IP Address.....Direct Connect Netmask:.....

Query/Retrieve supported:Y.....N

Storage Commitment :

SC IP address:.....SC Port Number:.....SC AE Title:.....

NETWORK HOST #3: Type:.....

Hostname:.....IP Address:.....Netmask:.....

Protocol:.....Port Number:.....AE Title:.....

Direct Connect IP Address.....Direct Connect Netmask:.....

Query/Retrieve supported:Y.....N

Storage Commitment :

SC IP address:.....SC Port Number:.....SC AE Title:.....

NETWORK HOST #4: Type:.....

Hostname:.....IP Address:.....Netmask:.....

Protocol:.....Port Number:.....AE Title:.....

Direct Connect IP Address.....Direct Connect Netmask:.....

Query/Retrieve supported:Y.....N

Storage Commitment :

SC IP address:.....SC Port Number:.....SC AE Title:.....

Installation and calibration

NETWORK HOST #5: Type:.....
 Hostname:.....IP Address:.....Netmask:.....
 Protocol:.....Port Number:.....AE Title:.....
 Direct Connect IP Address.....Direct Connect Netmask:.....
 Query/Retrieve supported:Y.....N
 Storage Commitment :
 SC IP address:.....SC Port Number:.....SC AE Title:.....

NETWORK HOST #6: Type:.....
 Hostname:.....IP Address:.....Netmask:.....
 Protocol:.....Port Number:.....AE Title:.....
 Direct Connect IP Address.....Direct Connect Netmask:.....
 Query/Retrieve supported:Y.....N
 Storage Commitment :
 SC IP address:.....SC Port Number:.....SC AE Title:.....

NETWORK HOST #7: Type:.....
 Hostname:.....IP Address:.....Netmask:.....
 Protocol:.....Port Number:.....AE Title:.....
 Direct Connect IP Address.....Direct Connect Netmask:.....
 Query/Retrieve supported:Y.....N
 Storage Commitment :
 SC IP address:.....SC Port Number:.....SC AE Title:.....

NETWORK HOST #8: Type:.....
 Hostname:.....IP Address:.....Netmask:.....
 Protocol:.....Port Number:.....AE Title:.....
 Direct Connect IP Address.....Direct Connect Netmask:.....
 Query/Retrieve supported:Y.....N
 Storage Commitment :
 SC IP address:.....SC Port Number:.....SC AE Title:.....

NETWORK PRINTER #1: Type:.....
 Hostname:.....IP Address:.....Netmask:.....
 Protocol:.....Port Number:.....AE Title:.....
 12 bits support:.....Life Size support:.....Color support:.....

NETWORK PRINTER #2: Type:.....
 Hostname:.....IP Address:.....Netmask:.....
 Protocol:.....Port Number:.....AE Title:.....
 12 bits support:.....Life Size support:.....Color support:.....

EC DECLARATION OF CONFORMITY

(Following the provisions of the medical devices directive 93/42/EEC)

We

Manufacturer

**GE Medical Systems SCS
283 rue de la Minière , BP 34
78533 BUC Cedex, France**

Manufacturing site

GE Medical Systems SCS
283, rue de la Minière, B.P. 34,
78533 Buc Cedex, France

Declare under our sole responsibility that the class IIa device :

Seno Advantage 2.2

Workstation, Workstation, diagnostic imaging, x-ray system

Ref. : addendum

GMDN Code: **40935**

Classification rule (93/42/EC Annex IX) : 10

To which this declaration relates is in conformity with the requirements of the medical devices directive 93/42/EEC which apply to it.

This conformity is based on the following elements :

- Information included in the documents :
Technical File Ref./ réf. : **DOC0732702** of the product to which this declaration relates.
- EC certificate: approval of full quality assurance system (annex II of the medical devices directive 93/42 EEC) delivered by G-MED (Notified Body 0459) on 21st of December, 2008 / Certificate N° 7840
- (If applicable) List of harmonized standards applied for CE marking
 - IEC60601-1-4 :2004
 - IEC60601-1-6 :2006
 - IEC60601-1-2 :2001
 - IEC60950-1 :2001

BUC, 22nd of March, 2010

Steve Slavens
Regulatory Affairs Director. MICT

This EC declaration of conformity supersedes the previous declaration dated November 3rd, 2009.

16

**ADDENDUM TO THE EC DECLARATION OF CONFORMITY dated 22nd of March, 2010**

Category	Identifier	Brief Description of Item
Basic System	5314329-x	Seno Advantage 2.2 HW SW
	5341143-x	Seno Advantage 2.2 - SW
	5340581-x	AW Volume Share4 core Software
	5363692-x	Seno Advantage right to use
	5183376-x	AW CONFIGURATION USB KEY
Hardware Options	5172850-x	19 Inches LCD Color monitor
	5310645-x	Disks Capacity Extension Option Package
	5344141 -x	RAM Capacity extension
	5183243 -x	HP 4GB Memory Package
	5161030-x	HP Workstation language Kit
	5135408 -x	EIZO 5 MP LCD MONITORS PACKAGE
	5165546-x, 5343105-x, 5343107 -x	EIZO 5 MP DIN6868-57 LCD MONITORS BUNDLE
	5257194-x	Barco Nio Monitor
	5318324 -x	Barco CORONIS Monitor
	5198838-x	Cable kit
	2224047-x	Key Pads for RWS
5183547-x, 5314326-x	Hardware options & spare parts	

ATTESTATION CE / EC CERTIFICATE

Approbation du Système Complet d'Assurance Qualité / *Approval full Quality Assurance System*
ANNEXE II point 3 Directive 93/42/CEE relative aux dispositifs médicaux
ANNEX II section 3 DIRECTIVE 93/42/EEC concerning medical devices

Fabricant (nom et adresse)
Manufacturer (name and address)

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

Mobile Senographe Essential

X-ray system, diagnostic, mammographic, stationary, digital

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
Declares under sole responsibility that the product **SecondLook® Digital (Product No: DSC002-01, DSC002-02, DSC002-03, DSC002-04, DSC002-05, DSC002-07, DSC002-08, DSC002-09, DSC002-10, DSC002-11, DSC002-12 and DSC002-19)** to which this declaration relates meets the essential health and safety requirements and is in conformance with the relevant EC directives listed below using the relevant section of the EC standards and other normative documents:

IEC 60950-1 (2001) First Edition
EN 55022:1994 +A1:1995 +A2:1997 Class A
EN 55024:1998 +A1:2002 +A2:2003
Medical Device Directive 93/42/EEC

MDD Class: Class I. The classification is based on the requirements of Rule 12 of Annex IX, of the Medical Device Directive.


Quality System Certification: ISO 13485:2003 as indicated on certificate number 9213, granted by Intertek.

CE Mark: Article 17 of the Medical Device Directive 93/42/EEC.

Date: January 26, 2010

iCAD, Inc., Authorized Representative in the European Community (as defined in Article 14 of the medical Device Directive: 93/42/EEC):

MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany
Phone: +49 511 62628630

Management Representative, John DeLucia
Vice President of Regulatory Affairs and Quality Assurance

Declaration of Conformity

Declared compliance: Council Directive EN 93/42/EEC
Standards UL 60601, CSA 60601, IEC 60601

Manufacturer: Barco N.V. MID
President Kennedypark 35
B-8500 Kortrijk
Belgium

Importer: GE Healthcare, Buc, France

Type of equipment: Medical Equipment, Class IIb

Model: Barco 5MP NON DIN LCD Grayscale Monitor

Serial number display 1: 1879045470

Serial number display 2: 1879045491

GEMS specification: 5257195 PSP rev 3

I, the undersigned, hereby declare that the equipment specified above conforms to the above declared Directives, standards and Requirements

Place: Kortrijk, Belgium

Date: 25/08/2010

Signature: 

Name: VERHAERKEN CHARLOTTE

Position: SYSTEM BUNDLING

**Declaration of Conformity
DBI™ Tables**

Hereby, Medical Positioning, Inc. declares that the #7402, 7407 and the #7425 DBI™ Tables are in compliance with the essential requirements and other relative provisions of Medical Device Directive 93 / 42 / EEC.

For evaluation of the compliance with these Directives, the following standards were applied:

- EN 60601-1
- EN 60601-1-2

The following importer / manufacturer is responsible for this declaration:

Company name: Medical Positioning, Inc.

Company Address: 1717 Washington Street, Kansas City, Missouri 64108, USA

Telephone: 816-474-1555

Facsimile: 816-474-7755

Web: medicalpositioninginc.com

Person responsible for making this declaration:

Michael G. Falbo, Sr.
Name (Full Name)

V.P. Product Development
Position / Title

December 28, 2005
Date

Legal Signature

EC DECLARATION OF CONFORMITY

Product: DBI TABLE, MODEL #'s 7402, 7407 & 7425

Manufacturer: MEDICAL POSITIONING, INC

Class: CLASS I, NON-MEASURING, NON-STERILE

DBI™TABLE 240 VAC Models

	DBI ^(d)	DBI ^(a)
Maximum Distributed Load:	226 Kg/ 400 lbs.	150 Kg/ 330 lbs.
Voltage:	240 VAC	230 VAC
Amperage:	2.6 A	1.0 A
Cycle:	50 Hz	50 Hz
Duty Cycle:	10%	10% (1 min/ 9 min off)

IEC 60601-1 Classifications:

- Class 1 Equipment
- Type B Applied Part
- Degree of Protection Against Ingress of Water / IPX0
- Equipment Not Suitable for Use in Flammable Anesthetic Mixture

All Electrical circuitry is isolated from chassis.

Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked Hospital Only or Hospital Grade

MEDICAL EQUIPMENT WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 2601-1 AND CAN/CSA C22.2 No. 601.1

Variants

Variants may include the above described surface incorporating ¹ Operation of the low voltage, electrically operated center column motor by means of a foot operated switch in place of the hand pendant. ² Addition of a second low voltage, electrically operated motor to enable Trendelenburg and Reverse Trendelenburg positioning of the surface and ³ A bearing supported intermediate section between the main section and the top of the low voltage, electrically operated center column motor to enable horizontal movement of the top surface while the base, center column motor, and casters remain in place.

Method of Manufacture

The product is manufactured in the following manner:

1. All steel components are cut, machined, welded, finish ground, primed and painted as sub assemblies.
2. All upholstered surfaces are prepared over structural composition board and fire retardant foam padding.
3. The steel components are final assembled, along with motor drives, casters and upholstered surfaces following the Assembly Procedures Document suitable for that product.

Product Operation And Use

The product provides a means of safe, comfortable, decubitus positioning of a patient for mammography procedures including, but not limited to, breast imaging and/or breast biopsy. One of the opening sections of the table is placed in the open position and the table is positioned so a portion of the mammography device can be positioned within the void created by the opened section. The portion of the mammography device positioned in the void is not attached to, or in contact with, the table surface. The table merely provides an opening large enough to accommodate the "buckey" of the device. The patient is positioned on their side, upon the surface of the table, facing the mammography device, with their breast positioned in a manner to be accessed by the mammography device. The remaining opening section of the table, on the side of the table away from the mammography device, is placed in the open position. The void created by this opening allow the physician and technologist greater access to the mammography device and aid in the correct positioning of the patient.

Risk Analysis Results


<u>Problem</u>	<u>Outcome</u>	<u>Risk</u>
No electrical power	Table can not be repositioned using motor drive(s)	None
Motor drive(s) not operational	Table can not be repositioned using motor drive(s)	None
Low voltage control not operational	Table can not be repositioned using motor drive(s)	None
Hand pendant not operational	Table can not be repositioned using motor drive(s)	None
Opening section can not be opened	Table can not be used	None
Opening section can not be closed	Latches must be repaired before using table	If table is used, patient could place hand in opening and possibly fall

Design Inspections


The product has successfully passed inspection and testing to the standards of UL 60601-1, IEC 60601-1 with Amendments No.1 and No.2 and CSA 22.2 No.601.1. by Underwriters Laboratories©

Labeling

The following is an accurate representation of the labeling to be affixed to the table. This label will be affixed to the center column motor.

Component Drawing

Design Drawing and Illustrations

Figure 1: A perspective view of a medical stretcher with a patient lying on it. The stretcher has a headrest and a footrest, and is mounted on a four-wheeled base.



Figure 2


Figure 5


Figure 4

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
Certification of Calibration

X-Rite, Incorporated certifies this instrument meets or exceeds published specifications applying to this instrument on the date it was calibrated. The instrument has been calibrated using standards whose accuracies are traceable to the standards listed in the products specification. The calibration is performed to these standards using equipment qualified at X-Rite to meet the specification for this instrument, and may use secondary instruments, and qualified illuminants.

Matrix 1

Measurement Delta as found:
 Luminance Ratio N/A
 Delta x N/A
 Delta y N/A

Measurement Delta as left:

Luminance Ratio 0.999
 Delta x -0.0010
 Delta y -0.0001

Matrix 3

Measurement Delta as found:
 Luminance Ratio N/A
 Delta x N/A
 Delta y N/A

Measurement Delta as left:

Luminance Ratio N/A
 Delta x N/A
 Delta y N/A

Matrix 2

Measurement Delta as found:
 Luminance Ratio N/A
 Delta x N/A
 Delta y N/A

Measurement Delta as left:

Luminance Ratio 0.999
 Delta x -0.0010
 Delta y -0.0001

Matrix 4

Measurement Delta as found:
 Luminance Ratio N/A
 Delta x N/A
 Delta y N/A

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