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Declaration of Conformity according to IVDD 98/79/EC

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Channel 3	37		
Channel 4	36,9		Y
Prewarming station	37		Y
10 Verify performance by running controls (ROTROL N) with reagents (mandatory):			
ROTROL N, lot number <u>41571001</u> , exp. date: <u>2014/03</u>			Values found accepted (y/n)
Gx [®] tem [®] (S), lot number <u>41706101</u> , exp. date: <u>2014/03</u>			Y
(star-tem [®] , lot number _____, exp. date: _____)			
First derivatives for all channels were checked for a smooth movement of the axis			Y
11 In case values of the performance run cannot be accepted:			
=> following corrective action was performed:			
12 Please check:			
X	X	Operating Manual has been read, especially Chapter 2, where important information is provided regarding risks, safety standards, operating conditions, and handling of reagents, samples and controls (mandatory)	

Final statement:

Mandatory checks are necessary to use the instrument for intended use. If one mandatory check cannot be marked please contact your service provider. After installation, the mentioned ROTEM[®] delta analyzer is in good condition and can be used for the intended use

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Please fax this form, the printout of the first derivative of the control measurement and the printout of the service values to Tem International GmbH at +49 89 46 42 95 22 or email a copy to support@tem-international.de.

	<ul style="list-style-type: none"> <input checked="" type="checkbox"/> The screen cable is connected on the rear side of the device to the VGA socket (mandatory) <input checked="" type="checkbox"/> The barcode scanner cable is connected to the screen (mandatory) <input checked="" type="checkbox"/> Touch screen power cable is connected with the power socket (mandatory) <input checked="" type="checkbox"/> USB interface at the rear side of the touch screen with the USB cable is connected with one of the USB ports of the ROTEM[®] delta (mandatory) <input checked="" type="checkbox"/> The keyboard is connected to one of the USB ports and placed on top of the device in front of the touch screen (mandatory) <input checked="" type="checkbox"/> The keyboard protection foil is placed on the keyboard <input checked="" type="checkbox"/> The country specific power cable is connected to the rear side of the ROTEM[®] delta and the power socket (mandatory) <input checked="" type="checkbox"/> The optional printer is connected to the USB port of the ROTEM[®] delta and to the power socket <input checked="" type="checkbox"/> The disc for cup ejector pin is placed around the cup ejector pin <input checked="" type="checkbox"/> The cup holders are placed at the temperature controlled area <input checked="" type="checkbox"/> The reagent rack is placed at the front of the working area <input checked="" type="checkbox"/> The accessory box is placed beside the ROTEM[®] delta <input checked="" type="checkbox"/> The ROTEM[®] delta is switched on (mandatory) <input checked="" type="checkbox"/> The printer (optional) is switched on <input checked="" type="checkbox"/> The touch screen is calibrated (Operation manual 6.7.1) <input checked="" type="checkbox"/> The software is in the right configuration (Operation manual 6.4) 															
5	<p>Take a cup&pin and place the pin on the shaft. Check if the pin is fitting too loose or too tight (mandatory)</p> <p>If the pin is fitting too loose => use the Axle Tool to spread the axis or use a small flat screw driver (very gently) If the pin is fitting much too tight => use a long (very gently) to press the axis together</p> <table border="1"> <thead> <tr> <th></th> <th>Pin fit (OK / loose / tight)</th> <th>Pin fit corrected (y/n)</th> </tr> </thead> <tbody> <tr> <td>Channel 1</td> <td>OK</td> <td>Y</td> </tr> <tr> <td>Channel 2</td> <td>OK</td> <td>Y</td> </tr> <tr> <td>Channel 3</td> <td>OK</td> <td>Y</td> </tr> <tr> <td>Channel 4</td> <td>OK</td> <td>Y</td> </tr> </tbody> </table>		Pin fit (OK / loose / tight)	Pin fit corrected (y/n)	Channel 1	OK	Y	Channel 2	OK	Y	Channel 3	OK	Y	Channel 4	OK	Y
	Pin fit (OK / loose / tight)	Pin fit corrected (y/n)														
Channel 1	OK	Y														
Channel 2	OK	Y														
Channel 3	OK	Y														
Channel 4	OK	Y														
6	<p>Check and note the service values in the service menu (mandatory):</p> <table border="1"> <thead> <tr> <th></th> <th>Variance found</th> <th>Variance accepted (y/n) (Target: < 12.000)</th> </tr> </thead> <tbody> <tr> <td>Channel 1</td> <td>2241</td> <td>Y</td> </tr> <tr> <td>Channel 2</td> <td>2021</td> <td>Y</td> </tr> <tr> <td>Channel 3</td> <td>2574</td> <td>Y</td> </tr> <tr> <td>Channel 4</td> <td>3582</td> <td>Y</td> </tr> </tbody> </table>		Variance found	Variance accepted (y/n) (Target: < 12.000)	Channel 1	2241	Y	Channel 2	2021	Y	Channel 3	2574	Y	Channel 4	3582	Y
	Variance found	Variance accepted (y/n) (Target: < 12.000)														
Channel 1	2241	Y														
Channel 2	2021	Y														
Channel 3	2574	Y														
Channel 4	3582	Y														
7	<table border="1"> <thead> <tr> <th></th> <th>Center found</th> <th>Center accepted (y/n) (Target: -832 to +708)</th> </tr> </thead> <tbody> <tr> <td>Channel 1</td> <td>17</td> <td>Y</td> </tr> <tr> <td>Channel 2</td> <td>-50</td> <td>Y</td> </tr> <tr> <td>Channel 3</td> <td>228</td> <td>Y</td> </tr> <tr> <td>Channel 4</td> <td>133</td> <td>Y</td> </tr> </tbody> </table>		Center found	Center accepted (y/n) (Target: -832 to +708)	Channel 1	17	Y	Channel 2	-50	Y	Channel 3	228	Y	Channel 4	133	Y
	Center found	Center accepted (y/n) (Target: -832 to +708)														
Channel 1	17	Y														
Channel 2	-50	Y														
Channel 3	228	Y														
Channel 4	133	Y														
8	<table border="1"> <thead> <tr> <th></th> <th>Amplitude found</th> <th>Amplitude accepted (y/n) (Target: 4402-4850)</th> </tr> </thead> <tbody> <tr> <td>Channel 1</td> <td>4642</td> <td>Y</td> </tr> <tr> <td>Channel 2</td> <td>4621</td> <td>Y</td> </tr> <tr> <td>Channel 3</td> <td>4665</td> <td>Y</td> </tr> <tr> <td>Channel 4</td> <td>4693</td> <td>Y</td> </tr> </tbody> </table>		Amplitude found	Amplitude accepted (y/n) (Target: 4402-4850)	Channel 1	4642	Y	Channel 2	4621	Y	Channel 3	4665	Y	Channel 4	4693	Y
	Amplitude found	Amplitude accepted (y/n) (Target: 4402-4850)														
Channel 1	4642	Y														
Channel 2	4621	Y														
Channel 3	4665	Y														
Channel 4	4693	Y														
9	<p>Check temperature (mandatory):</p> <table border="1"> <thead> <tr> <th></th> <th>Temperature found</th> <th>Temperature accepted (y/n) (Target: 35.5° - 38.5°C)</th> </tr> </thead> <tbody> <tr> <td>Channel 1</td> <td>36.9</td> <td>Y</td> </tr> <tr> <td>Channel 2</td> <td>37.0</td> <td>Y</td> </tr> </tbody> </table>		Temperature found	Temperature accepted (y/n) (Target: 35.5° - 38.5°C)	Channel 1	36.9	Y	Channel 2	37.0	Y						
	Temperature found	Temperature accepted (y/n) (Target: 35.5° - 38.5°C)														
Channel 1	36.9	Y														
Channel 2	37.0	Y														

Item	Quantity	Reference No
<input checked="" type="checkbox"/> eLine electronic pipette (mandatory)	1	700030
<input checked="" type="checkbox"/> eLine pipette filter	1	100048
<input checked="" type="checkbox"/> eLine pipette tips for the system pipette	1x96	400041
<input checked="" type="checkbox"/> Keyboard (mandatory)	1	400475-EN
<input checked="" type="checkbox"/> Keyboard protection foil	1	400476
<input checked="" type="checkbox"/> MC Rod pro (mandatory)	1	400046
<input checked="" type="checkbox"/> Set MC-Rod-Holder delta	1	200175
<input checked="" type="checkbox"/> Measuring cells (Cup&Pin pro) (mandatory)	1x20	200011
<input checked="" type="checkbox"/> Secure Viewer CD	1	200178
<input checked="" type="checkbox"/> Operating manual ROTEM [®] delta	1	200140-EN
<input checked="" type="checkbox"/> Pipette cable guide	1	500220
<input checked="" type="checkbox"/> Pipette and clean ring holder	1	200206
<input checked="" type="checkbox"/> Plastic clip for cup holder plus	8 (4 on holders)	400482
<input checked="" type="checkbox"/> Reagent rack ROTEM [®] delta	1	400484
<input checked="" type="checkbox"/> ROTEM [®] delta (mandatory)	1	200090
<input checked="" type="checkbox"/> ROTEM [®] delta power cable (mandatory)	1	700048-EN
<input checked="" type="checkbox"/> ROTEM [®] delta software	1	200140
<input checked="" type="checkbox"/> ROTEM [®] system USB Stick	1	400221
<input checked="" type="checkbox"/> Short User manual ROTEM [®] delta	1	200124-EN
<input checked="" type="checkbox"/> Short User manual ROTEM [®] delta SUR	1	200176-EN
<input checked="" type="checkbox"/> Touch screen Incl. power cable (mandatory)	1	400477
<input type="checkbox"/> Waste container	5x1	200133
<input checked="" type="checkbox"/> Parameter List	1	200086-EN
<input checked="" type="checkbox"/> Cleaning Ring Rod	1	400586
<input checked="" type="checkbox"/> Disc for cup ejector pin	1	400596
o		
o		
o		
o		

*If the printer is not delivered by Tem International GmbH, make sure to use a printer following IFC 60950 standards to comply with DIN EN 61010.

3 Please check:

- The on/off switch is always reachable (mandatory)
- The power cable can always be removed to disconnect the power socket in case of emergency
- ROTEM[®] delta is placed on an even, clean and robust surface/robust trolley
- The place where the equipment is set up is capable of supporting the weight of the system (22 kg) (mandatory)
- At least three earthed power supplies are available (mandatory)
- The equipment is not placed near open windows, radiators or other instruments that are producing heat or strong vibrations (mandatory)
- Sufficient ventilation of the ROTEM[®] delta is ensured, leaving a minimum distance of 20 cm between other instruments and walls (mandatory)
- The ventilation openings at the rear side of the equipment are not blocked (mandatory)
- Ambient temperature range of +15 °C to +30 °C (mandatory)
- Relative humidity range of 20% to 65% (mandatory)
- Installation site is in a height above sea level equal or lower than 2000 m (mandatory)

4 Please check:

- The system pipette is connected to the pipette connector on the right side of the ROTEM[®] delta (mandatory)
- The pipette and clean ring holder is connected to the right side of the device
- The pipette cable guide is pushed in the intended little borehole
- The touch screen is assembled
- The touch screen is placed on the ROTEM[®] delta

Installation Form for ROTEM[®] delta

tem[®]

Tem International GmbH advises to fill out the installation form for each instrument. The installation should be performed by trained personnel only.

Serial number of the instrument:	<u>2268</u>
Serial number of the monitor:	<u>L12.00 26620</u>
Serial number of the pipette:	<u>13002413</u>
IP Address:	_____
Gateway:	_____
Netmask:	_____
Gains:	_____
Software version:	<u>2.2.0</u>
Property of:	<u>ASI CRISTHO - Osp. S. Martino</u>

Printer Samsung CLP3
SN/Z4PRBAEBB006347

Installation at:	
Name:	_____
Contact person:	<u>Dressa Amelie Mulas</u>
Institute:	<u>ASI CRISTHO Osp. S. Martino</u>
Address:	<u>I.C.U.</u> <u>Via Carducci 35</u> <u>08170 CRISTHO</u>
Country:	<u>Italy</u>

Installation performed by:	
Name:	<u>Guido Bellussi</u>
Company:	<u>TEM International</u>
Date of installation:	<u>2013/05/2014</u>

1	Please check:																																
	<input checked="" type="checkbox"/> System consists of two packages (mandatory) <input checked="" type="checkbox"/> Number of additional packages for optional accessories <u>1</u> <input checked="" type="checkbox"/> No damage visible on the outside of the packages (mandatory) <input checked="" type="checkbox"/> Packages has been opened on the upper side <input checked="" type="checkbox"/> Cardboard with the accessories and the foam parts has been removed <input checked="" type="checkbox"/> Equipment has been taken out of the packages & placed on an even and clean surface <input checked="" type="checkbox"/> Accessories have been taken out of the packages (mandatory) <input checked="" type="checkbox"/> The pipette has been taken out of the white packaging (mandatory)																																
2	Check list of contents:																																
	<table border="1"> <thead> <tr> <th>Item</th> <th>Quantity</th> <th>Reference No</th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/> Accessories box</td> <td>1</td> <td>400472</td> </tr> <tr> <td><input checked="" type="checkbox"/> Adapter (<14 mm) sample preheater</td> <td>1</td> <td>400492</td> </tr> <tr> <td><input checked="" type="checkbox"/> Barcode scanner 3009 (mounted)</td> <td>1</td> <td>400566</td> </tr> <tr> <td><input checked="" type="checkbox"/> Barcode scanner 3009 holder (mounted)</td> <td>1</td> <td>400567</td> </tr> <tr> <td><input checked="" type="checkbox"/> Barcode scanner 3000 cable</td> <td>1</td> <td>400568</td> </tr> <tr> <td><input checked="" type="checkbox"/> Colour printer, USB (optional)</td> <td>1</td> <td>400222*</td> </tr> <tr> <td><input checked="" type="checkbox"/> Connecting cable eLine (mandatory)</td> <td>1</td> <td>400324</td> </tr> <tr> <td><input checked="" type="checkbox"/> Connection cable for printer, USB (optional)</td> <td>1</td> <td>400223*</td> </tr> <tr> <td><input checked="" type="checkbox"/> Cup holder pro (mandatory)</td> <td>4</td> <td>400016</td> </tr> </tbody> </table>	Item	Quantity	Reference No	<input checked="" type="checkbox"/> Accessories box	1	400472	<input checked="" type="checkbox"/> Adapter (<14 mm) sample preheater	1	400492	<input checked="" type="checkbox"/> Barcode scanner 3009 (mounted)	1	400566	<input checked="" type="checkbox"/> Barcode scanner 3009 holder (mounted)	1	400567	<input checked="" type="checkbox"/> Barcode scanner 3000 cable	1	400568	<input checked="" type="checkbox"/> Colour printer, USB (optional)	1	400222*	<input checked="" type="checkbox"/> Connecting cable eLine (mandatory)	1	400324	<input checked="" type="checkbox"/> Connection cable for printer, USB (optional)	1	400223*	<input checked="" type="checkbox"/> Cup holder pro (mandatory)	4	400016		
Item	Quantity	Reference No																															
<input checked="" type="checkbox"/> Accessories box	1	400472																															
<input checked="" type="checkbox"/> Adapter (<14 mm) sample preheater	1	400492																															
<input checked="" type="checkbox"/> Barcode scanner 3009 (mounted)	1	400566																															
<input checked="" type="checkbox"/> Barcode scanner 3009 holder (mounted)	1	400567																															
<input checked="" type="checkbox"/> Barcode scanner 3000 cable	1	400568																															
<input checked="" type="checkbox"/> Colour printer, USB (optional)	1	400222*																															
<input checked="" type="checkbox"/> Connecting cable eLine (mandatory)	1	400324																															
<input checked="" type="checkbox"/> Connection cable for printer, USB (optional)	1	400223*																															
<input checked="" type="checkbox"/> Cup holder pro (mandatory)	4	400016																															

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Il corso è stato tenuto da: Guido Bellussi (TECH INTERNATIONAL GmbH)